

BICO

TWO CHANNEL COAGULOMETER



- BICO is an high reliability coagulometer, able to perform the main coagulation tests.
- **PT**: result in Seconds, **RATIO**, **QUICK%**, **INR**. **aPTT**: two reagent test, with programmable incubation time, result in Seconds, **RATIO**. **TT**: result in Seconds and **RATIO**. **FACTOR**: result in Seconds, **RATIO** and **ACTIVITY%**. **FIBRINOGEN**: result in Seconds and concentration (mg/dl).
LCD display 40 columns X 2 rows. Membrane keyboard. 2 measurement channels. 44 sample incubation positions. 8 reagent positions. Serial interface RS232.
- Adjustable reading param. allows the use of any kind of reactive. 4 - 8 points calibration curves.
- Power supply 220V 50Hz 100VA Dimensions: 33 X 40 X 17h cm. Weight 8Kg.
Manufactured in according with $\text{C} \text{€}$ rules.



Dealer:



Manufacturer:



Minivolt Instruments S.r.l.

Via di Pietralata, 204/a

00158 Roma Italy

Tel. +39 06 418.20.89

Fax +39 06 450.42.75

www.minivolt.com

