


MONO

SINGLE CHANNEL COAGULOMETER



- MONO is an high reliability coagulometer, able to perform the main coagulation tests.
- 1 measurement channel. 2 reagent positions. 14 sample incubation positions.
- **PT**: result in Seconds, RATIO, QUICK%, INR. **APTT**: two reagent, with progr. Incubation time, result in Seconds, RATIO. **FIBRINOGEN**: result in Seconds and Concentration (mg/dl).
- LCD display 16 columns X 2 rows. Membrane keyboard. Serial interface RS232.
- Adjustable Reading param. allows the use of any kind of reactive. 4-8 points calibration curves.
- Power supply: 220V 50Hz 50VA Dim.: 11 X 20 X 11h cm. Weight 1.5 Kg.
- Manufactured in according with  rules.

Dealer:



Manufacturer:



Minivolt Instruments S.r.l.

Via di Pietralata, 204/a

00158 Roma Italy

Tel. +39 06 418.20.89

Fax +39 06 450.42.75

www.minivolt.com